


EXPORT ONLY



6.78MHz Monopolar RF
VOLNEWMER



VOLNEWMER®
6.78MHz Monopolar RF

- HIDDEN
EDGE TIP
- CURVED
SURFACE
- RADIAL
RF ENERGY
- CONTINUOUS
WATER
COOLING



T.+82-2-517-2114 | classys.com | info@classys.com
Images and texts are intellectual property of Classys. Copying of this material can be subject to charges of both civil and criminal law of legal justice. Copyright to Classys © All Rights Reserved



Volnewmer is the safest non-invasive monopolar radiofrequency device for skin tightening and lifting. It combines RF energy and cutting edge technology, giving you an overall skin treatment experience.

WHAT VOLNEWMER IS

COMBINING USER-FRIENDLY DESIGN AND TECHNOLOGY LEADS TO SIGNIFICANT IMPROVEMENT IN TREATMENT

VOLNEWMER is the safest RF treatment solution for skin lifting, tightening, and contouring by rejuvenation. It is the most advanced and innovative aesthetic system, with four interchangeable unique tip designed with a rounded surface tip to maximize RF power and treatment safety.

HIDDEN EDGE TIP

Maximized clinical results using side-effect suppression technology and safe Radiofrequency energy

CURVED SURFACE

Stable and most safe radial RF energy provided by the curved surface and rounded edge tip

RADIAL RF ENERGY

Radial RF energy delivers high level of treatment

CONTINUOUS WATER COOLING

Continuous water cooling prevents epidermal burns and able to deliver high level of RF energy




NEW 6.78MHz Monopolar RF

1

HIDDEN EDGE TIP DESIGN


The tip design is inward to prevent the phenomenon of energy clustering at the tip edges, ensuring safe energy transmission and minimizing side effects.




LESS PAIN AND MOST SAFE RADIAL RF ENERGY DELIVERED BY HIDDEN EDGE TIP

- The pain control enables higher treatment level with better clinical results.
- Ergonomics curved tips for perfect contact with the skin to optimize the RF energy delivery.
- Regular RF generating on the whole surface avoids adverse effects on the epidermis.


* The RF energy remains constant at the same level for up to 60 seconds.




After 15 sec



After 30 sec




After 45 sec



After 60 sec

VOLNEWMER®



2

CURVED SURFACE & RADIAL RF ENERGY

RF energy delivers uniform energy to the dermis by ergonomics curved tips. The curved surface of the tip makes perfect contact between the tip and the skin. The RF thermal energy of VOLNEWMER forms radially, enabling the treatment at a high level.



COLLAGEN RE-STRUCTURE AND RE-GENERATION BY OVERLAPPING ENERGY

- RF energy renew the collagen, promote skin barrier, improve elasticity and tighten the skin.
- Accelerates collagen regeneration and then enhances ECM.

3


CONTINUOUS WATER COOLING SYSTEM

VOLNEWMER employs a stable high-frequency output technology to boost energy delivery efficiency and reduce pain and side effects by emitting a single energy pulse per REP, instead of multiple pulsations. The device also includes an optimized cooling system designed for efficient RF energy production, working through continuous water cooling as skin resistance increases.



PROVIDING QUICK AND COMFORTABLE TREATMENT

- Continuous water cooling alleviates pain with excellent safety and directly protect epidermis injury during RF generation.



PROVIDING STIMULATES COLLAGEN REGENERATION IN THE PAPILLARY DERMIS

- Hydra water cooling affects the heat energy formation in the papillary dermis, and it produces immediate effects right after treatment.





EASIER FOR PHYSICIAN WITH BETTER CLINICAL RESULT

1

OPTIMAL REPEAT MODE

- Repeat mode enables faster treatment.
- VOLNEWMER's hidden edge and curved tip surface generate the radial flow of the electrodes and create overlapped thermal energy formation for optimizing repetitive, overlapping procedures.



2

2 HAND-PIECES

It enables the application of quick and simple treatment. The two handpieces reduce the hassle of changing tips during treatment.





NEW 6.78MHz Monopolar RF

3

FLEXIBLE & EASILY TILTING TIP

The tip head of VOLNEWMER is designed to tilt in all directions, allowing it to move according to the curvature of the treatment area, thereby enhancing the stability of contact between the tip and the skin.



4

REAL TIME IMPEDENACE DETECTING SYSTEM

Volnewmer detects the change in impedance value across the entire tip during the shot and contains a technology that stops the RF generating if even a part of the tip detached the contact surface to prevent side effects.

VOLNEWMER®

5

NO TIME LOCK

The Treatment schedule should be determined by physician's judgment. No time limit on consumables means physicians can combine up to four types of tips during one procedure.


6

4 TYPES OF TIPS TO TREAT MULTIPLE AREAS

Volnewmer's different types of tips are designed to allow customizable treatment for face and body.



Various procedural designs are available such as face (eye & forehead), Body and Décolleté



ARMS

ABDOMEN

KNEES & THIGH

BUTTOCKS


NEW 6.78MHz Monopolar RF

CLINICAL RESULTS

FRONT VIEW



BEFORE TREATMENT




8 WEEKS AFTER




12 WEEKS AFTER


VOLUME CHANGE



BEFORE TREATMENT




8 WEEKS AFTER




12 WEEKS AFTER

* Volume change visualizes quantitative analysis of volume changes (a total of 17.43 cc volume increase 3 month post treatment)


VECTOR CHANGE



BEFORE TREATMENT



8 WEEKS AFTER




12 WEEKS AFTER

* Vector change visualizes the degree of skin lifting with the arrow

VOLNEWMER®


CLINICAL RESULTS

VOLUME EFFECT




BEFORE TREATMENT

8 WEEKS AFTER



* Volume change visualizes quantitative analysis of volume changes



* Vector change visualizes the degree of skin lifting with the arrow


LIFTING EFFECT



8 WEEKS AFTER

12 WEEKS AFTER

LIFTING & WRINKLE REDUCTION



BEFORE

12 WEEKS AFTER

BEFORE

12 WEEKS AFTER

NEW 6.78MHz Monopolar RF


CLINICAL RESULTS

SKIN ANALYSIS



BEFORE TREATMENT

12 WEEKS AFTER




SKIN ANALYSIS



BEFORE TREATMENT


12 WEEKS AFTER



* Higher percentage indicates more skin improvement on each category

VECTRA 3D CAMERA REPORT

PERCENTAGES OF PEOPLE WHO IMPROVED THE SKIN BY EACH CATEGORY



VOLNEWMER®